

Microscopy and microscopes

General Microbiology - Laboratory

Cañada College - Fall 2008

Instructor: Tamas Torok, Ph.D.



Today's agenda

- **Bright-field light microscopy**
 - history of microscopy
 - today's light microscopes
 - structure and function of the bright-field light microscope
- **Other microscopes**

Ernest Abbe (1840-1905)

- **“... a non-self-luminous particle, which is illuminated by an extraneous source, gives rise to diffracted light rays...”**
- **“...to form a good microscopic image as many of the diffracted rays as possible should be intercepted by the objective...”**

History of the light microscopy



Leeuwenhoek's microscope





Fig. 1 - Microscope made by Anton van Leeuwenhoek in VIIth century.

State-of-the-art



Abbe's formula

- Resolving power (d)
- Wavelength of light (λ)




Figure 2

- Index of refraction (n)
- One half the angular aperture (μ)
- Numerical Aperture
- $NA = n (\sin \mu)$
- $d = \lambda/NA = \lambda/n (\sin \mu)$





Figure 1

Numerical aperture and resolution

- In confocal and fluorescence microscopy, the resolution may exceed these limits
- Other factors, such as low specimen contrast and improper illumination may lower resolution



Objective information

60x Plan Apochromat Objective




Figure 1

	NA	WD [mm]
10x	0.45	4.00
40x	0.95	0.14
100x oil	1.40	0.13

NA	.25	.50	.95
Depth [μm]	8.0	2.0	.10

Enhancing resolution power

Oil Immersion and Numerical Aperture





Figure 1

- Numerical aperture is related to the imaging medium
$$NA = n(\sin \mu)$$
- Because μ cannot be greater than 90° ($\sin \mu = 1$), the maximum possible numerical aperture is determined by the refractive index (n) of the immersion medium, including water ($n = 1.33$), glycerin ($n = 1.47$), immersion oils ($n = 1.52$)
- $d_{th} = 0.25 \mu m$ ($\lambda = 550 \text{ nm}$; $NA = 1.40$)

Other light microscopic techniques

- **Differential interference contrast (DIC)**
 - a mode of contrast generation in microscopy that uses polarized light and an “analyzing filter”, yielding an image with a shadow relief
- **Dark field**
 - illuminates the specimen but does not admit light directly to the objective
- **Stereo dissection**
 - perceived depth by transmitting twin images that are inclined by a small angle to yield a true stereoscopic effect

Phase contrast microscope



- Phase of the light diffracted by the specimen is altered by approximately $1/4$ wavelength
- Human eye is sensitive only to color (light frequency) or to light intensity (wave amplitude)
- Rings cause (destructive or constructive) interference - better contrast
- Zernike received Nobel Prize in 1953

Fluorescence microscope

- **Excitation light irradiates the specimen**
- **Fluorescent light is separated from the brighter excitation light and only the emission light reaches the eye**
- **Fluorochromes are specific in their attachment to biological structures**
- **Energy-transfer fluorochromes**





Figure 3


Laser scanning confocal microscope




- Resolution: $0.2\ \mu\text{m}$
- Optical sectioning results in 3D image
- Limited by scanning speed



X-ray microscope




Transmission electron microscope



- Electron wave ($\lambda < 0.005 \text{ nm}$)
- Thousand fold increase in resolution, hundred fold in depth of field
- High vacuum needed
- Extensive sample preparation
- Sample is always dead, images often contain artifacts


Scanning electron microscope

- Magnification ranges from 15x to 200,000x
- Resolution: 5 nm



©Eric B. Workman

Atomic force microscopy



- Image at atomic resolution
- In the non-contact mode (distance $>10\text{\AA}$ between the tip and the sample surface), Van der Waals, electrostatic, magnetic or capillary forces produce images of topography
- In the contact mode, ionic repulsion forces take the leading role

A Window on the Future

Consider gallium nitride, a semiconductor whose promising applications include light-emitting diodes that produce a rainbow of colors: until recently it was impossible to make atomic-resolution images of such a material. Under a transmission electron microscope, heavy atoms like gallium scatter electrons so much that they swamp the signal from lightweight neighbors like nitrogen. Then the One-Angstrom Microscope (OAM) debuted at Berkeley Lab's National Center for Electron Microscopy (NCEM).



Modern Microscope Component Configuration




Figure 1