

Verification of the Bubble Column Method

The bubble column method was tested by comparing the K value obtained from the bubble column method through Eq. (5), with the K value from surface tension measurements through Eq. (6), for dilute surfactant solutions of sodium dodecyl benzene sulfate (SDBS).

Rearranging Eq. (4) gives

$$K = \frac{AD}{afz} \ln \frac{C(z)}{C_b} \quad (5)$$

Combining Eq. (1) and Eq. (3) gives

$$K = \frac{-1}{RT} \frac{d\gamma}{dC} \quad (6)$$





Fig. 3. Testing the bubble column method on SDBS solutions. (a) $K = 56 \pm 3 \mu\text{m}$, from surface tension measurements. (b) $K = 56 \pm 13 \mu\text{m}$, from the bubble column method. This agreement validated the bubble column method.